

CaCl₂ Chemical Competence (for *E. coli*)

DAY I: STREAK OUT

DAY II: O/N CULTURE

1. Inoculate 2 mL O/N culture from single colony from part I

DAY III: GROW

2. Subculture (back-dilute) 1:25-1:50 in 100mL fresh media
2mL culture □ 100 mL media for 102 mL total is fine
3. Grow cells to mid-log (OD₆₀₀=0.3-0.5)
Put in shaker at 200 rpm at 37C
Grow for ca. 1.5 hours
4. Harvest cells by centrifugation, pour off supernatant
4000 rpm for 20 min, or 5000 rpm for 10 min; more time as needed
Use two 50 mL Falcon conical tubes
Fast-cool the centrifuge to 4C
5. Resuspend cells in COLD (4C), STERILE CaCl₂
Put 10 mL CaCl₂ in each conical, for a total of 20 mL
Use 100 mM (.1 M) CaCl₂ for most cases
Use 50 mM (.05 M) CaCl₂ for M13 hosts, i.e. TG1 or JM103)
6. Incubate O/N in the cold room

DAY IV: FINISH

7. Pellet cells by centrifugation, pour off supernatant
4000 rpm for 20 min, or 5000 rpm for 10 min; more time as needed
Fast-cool the centrifuge to 4C
8. Resuspend cell pellet in COLD (4C), STERILE CaCl₂
Put 2 mL CaCl₂ in each conical, for a total of 4 mL
Use 100 mM or 50 mM as needed

If using immediately: store on ice for 2-12 hours prior to transformation; after 24 hours, efficiency of transformation drops.

If storing: add 0.5 mL COLD (4C), STERILE 50% glycerol to each tube, for 1 mL total. Aliquot 310 mL into eppi. tubes and store at -80C.