

Polyacrylamide Gel Electrophoresis

- (1) Resuspend and lyse the cell pellets from the induction.
 - Resuspend the pellet in 50 mL TE
 - Add 50 mL 2X Sample Buffer (SB)
- (2) Temperature shock.
 - Alternate 5 min. heating in the 100C block and freezing in the -80C freezer.
 - Heat, freeze, heat, freeze, heat. (Beware! Lid can pop open during heating.)
 - End on heat, because freezing can cause unrelated proteins to “globulate.”
 - Let cool to R.T. before loading.
- (3) Load the solution.
 - The sample with the lowest OD should be the 10mL (max. vol.) sample
 - Use vol. and OD in $MV=MV$ to standardize the amount of protein being loaded.
 - Re-freeze unused solution at -80C.
- (4) Run at 100V until the dye hits the separating gel, then crank up to 160V.
- (5) Stain.
 - Remove, put in Gladware with ventilated lid.
 - Cover in Fairbanks stain (contains Coomassie Blue).
 - Microwave for 40-50s.
 - Put Gladware on shifty thing in hood for 10 min.
- (6) Destain.
 - Pour stain in waste bin.
 - Add Fairbanks destain.
 - Microwave for 50s.
 - Add two small Kimwipes.
 - Put Gladware on shifty thing in hood.
 - Check after 30 min-1 hr to preview the gel.
 - Remove Kimwipes, repeat first five destain steps, let destain overnight.

Setting up the gel rig:

	<ol style="list-style-type: none">1. Lock gels into place. Make tops flush with bottom ridge on green padding.2. Mix buffer soln: 450 mL H₂O + 50 mL 10X
--	--

3. Fill inner well to top first. Push off bubbles
Wait to see if it leaks.

4. Pour remainder of soln into outer well.