


## Using the Spectrofluorimeter

### POWERING ON AND SETTING UP:

1. Turn on cooling box (1)  
Let run for 5-10 minutes
2. Turn on incubator (2) switches in the following order:  
Lowest switch on bath; highest switch on box; on/off button on box  
Let incubator run for ~5 minutes
3. Turn on computer (8)  
\*Do not open Olis Globalworks program until AFTER photon counter and control box are on
4. Turn on "Power" to control box (3)
5. Push and hold "Ignite" on control box; hold until green light (4) comes on  
You want the control box to display 147-148 watts steadily
6. Turn on photon counter (5)  
\*\*\* THE COUNTER SHOULD NOT BE EXPOSED TO AMBIENT LIGHT  
Make sure all four manual slits, ESPECIALLY the one near the photon counter, are pushed in so no ambient light can enter and damage it.
7. Turn on Olis computer (7)
8. Place the cuvette holder in the well (6)



1. Cooling box	5. Photon counter
2. Water bath	6. Well
3. Control box	7. Olis computer
4. Green light	8. Computer computer

# Using the Spectrofluorimeter

## DATA COLLECTION SETUP

*Program:* Olis Globalworks

Click Data Collection tab

Open “DM45”

If: “The masterboard has timed out,” then: “Retry” until functional.

*Menus:*

1) Parameters

Where you change hardware settings

Two PMTs and the photon counter (PC)

“Hardware reset and calibration” if stuff isn’t working

1. Operational Modes

Where you can tell it about the assay

Data Reduction Mode: Always “Photon Counter”

Data Collection Mode: “Assay” (time-based)

2. Live Display

Set “emission” = “excitation” = 310 nm

Set duration of assay to 300s for FtsZ

Set number of data points and frequency of data point collection, usually every 0.25 s

\*make sure (# data points)(frequency data taken) = (duration of assay)

*Saving Files:*

At the end of a reading, the Globalworks should automatically prompt the user to post data.

If it does not, click “Post these data” in the top right of the program window

Select the trace by clicking on it

Right-click the desired trace and save as an ASCII file **Using the Spectrofluorimeter**

## DATA COLLECTION

*Loading and Reading*

1. Pre-chill all buffers
2. Using a gel loading tip, wash out the cuvette well with 2-3mL of water before initial use; take care to fill the cuvette fully and not to

- overflow it. Hold the cuvette by the top and bottom, or by the corners.
3. Wash out the cuvette with buffer using gel loading tips 2X between runs.
  4. Load all buffers, load protein, pipet well to mix.
  5. Load any other proteins, pipet well to mix.
  6. Final volume must not exceed 200  $\mu$ L; this includes later additives like GTP
  7. Wipe off cuvette with lens paper
  8. Load into spectrofluorimeter, with "Q" on cuvette facing the left (for consistency)
  9. On the computer, click "Collect Data," then press the spacebar
  10. Let the first 30-60s pass to get a baseline reading; in the meantime, empty the old cuvette
  11. After 30-60s, add 2 $\mu$ L GTP directly to the cuvette in the spec and mix well  
Brad uses a gel loading tip; Dan and Amy switch to a larger tip to mix  
Do not put the tip too far into the cuvette as this will affect the reading
  12. Wash the old cuvette during the remainder of the assay