

## DpnI mediated site-directed Mutagenesis

A highly effective simple method for making site directed lesions in plasmids without subcloning based on the work of Fisher and Pei (1997).

### **1. Amplification of mutant DNA**

DNA template plasmid 5-20 ng  
10x pfu DNA polymerase buffer 5.0  $\mu$ l  
25uM oligo 1 0.5  $\mu$ l  
25uM oligo 2 0.5  $\mu$ l  
10mM dNTP 1.0  $\mu$ l  
Pfu DNA polymerase (2.5 units) 1.0  $\mu$ l  
fill w/ddH<sub>2</sub>O to 50  $\mu$ l

#### PCR conditions

95° 30 seconds,  
18 cycles of : 95° 30 seconds, 55° 1 minute, 68° 2 minutes/kb of plasmid length

### **2. Degradation of methylated (parental) DNA with Dpn I**

Cool down PCR reaction.  
Add 1 $\mu$ l Dpn I (10 unit) to PCR reaction 37°C and incubate for 1 hr.

### **3. Transformation into *E. coli***

Place 200  $\mu$ l competent cells (1 x 10<sup>8</sup>/ug efficiency or greater) DH5 $\alpha$  on ice.  
Add 1-2  $\mu$ l of digested PCR reaction.  
Incubate on ice 15-20 minutes.  
Heat shock 45 seconds, and return to ice for 2 minutes.  
Add 1 ml LB broth. Incubate at 37°C for 1 hr.  
Plate 100  $\mu$ l aliquot on appropriate antibiotic plate.

### **4. Miniprep 6 colonies and digest plasmids looking for mutant.**

### **5. Sequence the lesion to confirm the change.**

#### **Notes:**

##### **Design of oligos**

Oligos should be perfectly complementary (actually haven't checked out if this is required). The oligos should contain 12-15 bp on each site of the lesion. It is useful to engineer a restriction site addition or loss into the oligos for ease of screening. We have used the protocol for deletions using 17 bp on each side of the deletion.

##### **Plasmid size**

We have used the protocol for plasmids over 18 kb in length. Plasmids below 10 kb seem to work routinely. In most cases, a vast majority 80-100% of colonies are correct. For longer plasmids it may be necessary to work with PCR to optimize synthesis.

##### **Introduction of lesions**

Note that this is a linear amplification, so that the template only gets copied one time. Hence the odds of introducing lesions when using a high fidelity polymerase are very low.

##### **Analysis of PCR product by gel electrophoresis**

Run an agarose gel ( 5 -10  $\mu$ l). If the product is visible, the mutagenesis will almost certainly work. Even if the product is not visible, one often still gets the mutant.

### **Using other polymerases**

We have started using TaKaRa Ex-Taq instead of Pfu polymerase with great success using large plasmids (15-20 kb size). While getting the amplification to work efficiently with Pfu take a bit of fiddling, amplification with Ex-Taq appears very robust for large plasmids.

### **References**

Fisher, C. L., and Pei, G. K. (1997). Modification of a PCR-based site-directed mutagenesis method. *Biotechniques* 23, 570-574.