

## A PCR approach to characterizing MosSCI

Make high quality genomic DNA. Our protocol for doing this can be found at <http://neurosci.wustl.edu/nonetlab/resources>.

We use PCR to amplify up inserts

for II

NM3880 5' AGGCAGAATGTGAACAAGACTCG

NM3887 5' ACCGGAAACCAAAGGACGAGAG

NM3888 5' ACGCCCAGGAGAACACGTTAG

NM3884 5'ATCGGGAGGCGAACCTAACTG


for IV

NM3881 5'CAAACGGAGCACCAGGAAAGC

NM3885 5' AAACTCCAAACACACCAGTCAC

NM3889 5' CCAAACAAGTGTCGTTGACCCAG

NM3890 5' CATATCCGCCAAGGACGCTC



PCR reactions are performed using NEB LongAmp polymerase under the following conditions

#### Reactions

5 ul 5 X LongAmp Buffer

0.75 ul 10mM dNTPs

0.5 ul 3884 oligo 10 uM

0.5 ul 3880 oligo 10 uM

5 ul template ( of a 1:50 dilution of our Genomic DNA prep ~ 100 ng total)

Water to 25 ul

1 ul LongAmp polymerase.

#### PCR conditions

30 second 95°C

35 cycles of

94°C for 10 seconds

60°C for 50 seconds

65°C for 12 minutes +10 second per cycle

65 for 10 minutes

We use an MJ Research PTC100 thermocycler using the EXTEND function to increase the extension time each cycle. I doubt this is really needed.

#### Samples

	PCR size bp	
	3880/3884	3880/3888
1 N2	2942	1674
2 EG4322	4235	2967
3 EG4887 myo2 cherryH2B unc543'	8453	7185
4) NM2237 plasmid multicopy insertion	-----	6105
5) NM2237 plasmid single copy insertion	7373	6105
6) NM2237 plasmid single copy insertion #2	7373	6105



- 3) The lane 7 and lane 8 integrants actually look indistinguishable under the microscope looking at GFP fluorescence, but it is likely that the #7 has some unusual structure. This approach can distinguish such irregularities.
- 4) We also have PCRed using oligos in vector sequences and shown that the multicopy insertion amplifies with oligos in the Amp and Ori, while all the other strains do not.
- 5) One could also look for the presence of Mos1 transposons that may have hopped back in using oligos. This is short range PCR and does not require a high template quality.
- 6) We have not succeeded in PCRing up from single worm lysis (or multi-worm lysis) for these long range PCR reactions.
- 7) We have not attempted to optimize the oligonucleotides for long range PCR. Barnes (1994 PNAS 91, 2216) claims that longer oligonucleotides that anneal at 68°C work more robustly for long range PCR. For LongAmp taq which recommends 65°C extensions, oligos which can permit the melding of the anneal and extend steps may improve performance further.
- 8) We have tried the following polymerases: ExTaq (works up to ~ 10kB well) using 60°C anneal, 72°C 10 minute extensions. KlenTaq LA is less effective than ExTaq using their pH9.0 buffer and recommended annealing and extension times.
- 9) You may think of using LongAmp for other purposes since it gives such robust products. Beware we have found the error rate is extremely high!!
- 10) If you don't get a product using either approach (across the whole insert, or from one arm to outside, one can still confirm there is an integration using one one oligo outside the arm and an oligo in unc-119 (one one side) and one oligo outside the arm and one oligo in your gene for the other side. These PCRs can be designed to create small product (<2 Kb) and done on single worms before high quality genomic DNA is even prepared.